

AD619496

OPTICAL PERFORMANCE
of
MULTI-CHANNEL UV and VISIBLE SPECTROGRAPH

COPY _____	OF _____	<i>2</i>
HARD COPY	\$. 2 . 0 0 -	
MICROFICHE	\$. 0 . 5 0	

by

JAMES J. DEVLIN, S. J.

21P


Boston College, Chestnut Hill, Massachusetts, 02167

CONTRACT NO. AF19(428)-4305

Project No. 6327

FINAL REPORT

27 May 1965



Period Covered 1 April 1965 to 30 June 1965

Prepared for

Air Force Cambridge Research Laboratories
Office of Aerospace Research
United States Air Force
Bedford, Massachusetts

This research was supported by the Advanced Research Projects Agency under
ARP^A Order 254

ARCHIVE COPY

**BEST
AVAILABLE COPY**

Requests for additional copies by agencies of the Department of Defense, their contractors, or other government agencies should be directed to:

Defense Documentation Center (DDC)
Cameron Station
Alexandria, Virginia 22314

Department of Defense contractors must be established for DDC services or have their "need-to-know" certified by the cognizant military agency of their project or contract.

All other persons and organizations should apply to the:

Clearinghouse for Federal Scientific
and Technical Information (CFSTI)
Sills Building
5285 Port Royal Road
Springfield, Virginia 22151

INTRODUCTION

The spectrograph under test is an Ebert Mount Spectrograph designed and built by ITT Research Institute, Technology Center, Chicago, Illinois.

The sources used in testing this spectrograph were Helium and Neon Geisler tubes, powered by a 3,000-Volt transformer, a low pressure mercury tube with a quartz wall and a Mazda tungsten filament lamp. The Mercury sources were powered by a 110-Volt A.C. power source, and the Mazda lamp by a 6-Volt transformer.

Throughout the tests a 103 Fa emulsion on 35mm. film was used. This emulsion was chosen because it is highly sensitive in the visible range and reasonably sensitive in the longer wavelength range of the UV.

A standard development procedure was followed in all tests. The film was developed in D-19 at 20° C for five minutes and fixed for ten minutes.

SPECTRAL RANGE

A Hg lamp source registered the Hg line 2536 Å at the lower end of the spectrogram. A Neon Geisler tube registered the line 6382 Å at the upper end of the spectrogram. An extrapolation showed that the spectrogram would record lines in the wavelength range from 2500 Å to 6500 Å.

PLATE FACTOR

A number of line pairs in the Mercury and Neon spectrum were identified and the distance between these measured. The average plate factor was determined to be 111 Å/mm .


RESOLVING POWER

By controlling the exposure time of the Hg source it was possible to identify the Hg lines 3023.5 Å and 3021.5 Å . These lines were clearly resolved, verifying the resolving power as 2 Å .

FIELDS OF VIEW

The field of view was determined by mounting a He tube on a goniometer and setting the tube at a series of angles with reference to the optical axis of the fibers. No optics were used between the He tube and the fiber ends. The Helium tube has a very narrow bore thus providing a source width of approximately a millimeter. The exposure time was standardized at 20 seconds for those tests.

The field of view is defined by the angle θ subtended at the fiber axis by the source.



Two lines in the Helium source were selected for study 5875Å and 3888Å. The density of these lines at successive angle positions of the source were measured. The relative intensities corresponding to these densities were determined from the characteristic curves and plotted against the angles on semilog paper. The lack of symmetry in these curves can be ascribed to the fact that it is extremely difficult to align the fiber axis on a fixed position. The field of view of the instrument could have been determined without the fibers but this would not be representative of the instrument as designed for flight operation. The data indicates a field of view with relatively high intensity for about 6° with a lower response out to 12° . The data is plotted in figures I through IV.

SPEED

The speed of a spectrograph is dependent on many factors. In addition the term is not well defined.

The most reliable source for testing the spectrograph would be a black body operating at a temperature of 3100K or 2900K. Since such a source was not available at this laboratory an automobile headlamp bulb with a tightly coiled tungsten filament and an input of 16 watts at 6 volts was used. The diameter of the filament was 0.48mm. and its length 5mm. In terms of the conventional definition of brightness the energy output is 6.66 watts/mm^2 .

This source within a reasonable error is presumed to be a black body radiator. An exposure of one second was taken with the Mazda lamps

located at a distance of 40 and 30 cm respectively. The film was developed under the standard conditions given above. A microphotometer trace was taken of the transmission T vs wavelength, for the exposures. The transmission values were converted to densities, $D = \log \frac{100}{\%T}$. The results are plotted in figures V and VI.

The characteristic curves of this emulsion were obtained by exposing the 103 Fa emulsion to the radiation from a D. C. arc operating at 5 amps and using 3/8 in. iron electrodes. The ITT spectrograph does not have sufficient dispersion to record the lines of the iron arc. This exposure was made on a Jarrell Ash 1.5 meter Wardsworth mount spectrograph. The DC arc radiation was attenuated at the slit by a rotating sector with seven steps in the ratio of 1.585 to 1. The exposure time was of the order of 5 seconds so as to produce a reasonable level of densities.

By the usual two-step method the data was reduced to provide a series of curves relating the optical density to the log of the intensity I . This relation is a function of the emulsion and not of the spectrograph. The intensities given on these curves are relative.

Included in this report is a graph of the black body radiation for a source at 2900°K taken from the Tables of Black Body Radiation Functions by Mark Pivovorsky and Max R. Nagel. The ordinate is given in watts per cm^3 per steradian.

A fiber bundle consists of four fibers of 0.01 in=0.254mm. The intensity per steradian picked up by a single fiber would depend on the position of the source. The solid angle could then be calculated. Since the exposure $E = \int I dt$ the total energy recorded would then be calculated. This calculation where, I = intensity and t = time, would be reliable in the case of a single line recorded on the spectrogram. If a continuum is recorded, a correction must be made for slit width. The procedure may be found by referring to Chapter Six of PRACTICAL SPECTROSCOPY by Harrison Lord and Loofbourov, or some other standard reference book.

The abscissa is in milimicrons. Thus by choosing a band width the intensity in watts per cm^2 per steradian may be determined for any desired wavelength in the spectrum. The density that corresponds to this wavelength may be determined from the density vs wavelength curve. Then the relative intensity on the characteristic curve may be converted to an absolute value.

There is a source of error always present in this instrument, namely, the alignment of the fiber ends to the source. This is unavoidable since the instrument was designed to bend the fibers and direct them towards several sources simultaneously. However, the relatively wide field of view allows for the recording at a spectrum of a source when the fiber terminal is pointed in the direction of the source.


FIBERS


Several exposures were run off with each fiber bent at various angles within the limit of safety of the fibers. A Helium Geisler tube was used for these tests. The spectra compared favorably with those used to determine the angle of view.


RECOMMENDATIONS


No mechanical tests were run to determine the ruggedness of the instrument. However, the open reels should be replaced by light proof magazines.


Whenever this spectrograph is used one fiber should be directed toward a quartz tube mercury lamp of moderate brightness in order to record a calibration spectrum for wavelength determinations.















































RELATIVE INTENSITY

FIG. III FOR $\lambda = 5875 \text{ \AA}$

